




- ◆ Dual channel
- ◆ 3 Modes: Standard, Modulate and Burst
- ◆ Pulse Width: Adjustable: 40 - 260 Microseconds
- ◆ Timer: 60 minutes, 30 minutes and continuous
- ◆ Pulse Amplitude: Adjustable, 0-80mA each channel or -20mA peak maximum negative
- ◆ Wave Form: Asymmetrical Bi-Phasic square pulse
- ◆ System Includes: **Stimulator, Lead Wires (New FDA compliant), 4 Self-stick Hypoallergenic Reuseable Electrodes, 9 Volt Battery, Hard Plastic Carrying case, Instruction Booklet.**

ProM-300 CONTROLS AND INDICATORS



TECHNICAL SPECIFICATION

ProM-300 (TENS THREE MODS WITH TIMER)

Channels:	Dual, isolated between channels
Wave form:	Modified square with zero net direct current (DC) component
Pulse Amplitude:	Constant current 0 to 80 mA each channel, adjustable (500 ohm)
Pulse Frequency:	2 Hz to 50 Hz Adjustable
Modulation Mode:	Pulse width us automatically varied in a cycle pattern over an interval of nominally 4.0 second
Burst Mode:	9 +/- pulses at maximum pulse rate twice every second
Maximum charge:	20 micro coulombs per pulse
Power Source:	9 v alkaline battery or similar rechargeable sell
Size:	23 x 64 x 91 mm
Weight:	130 grams (including battery)

All values have 20% +/- tolerance

Voice: 215.938.0200 ~ Toll-Free 1.800.472.9190 ~ Fax: 215.938-9309

www.ProMedSpecialties.com ~ 1800 Byberry Rd. Unit 905 ~ Huntingdon Valley, PA 19006 ~ USA